

ERA

Expires February 28, 2017

Analyte	Method		Analyte	Method
Metals			Synthetic Organic Chemicals	
Antimony	200.8		Benzo(A)Pyrene	525.2
Arsenic	200.8		Carbofuran	531.1
Barium	200.8	200.7	DBCP (1,2 Dibromo-3-Chloropropane)	504.1
Beryllium	200.8	200.7	Di(2-Ethylhexyl)Adipate	525.2
Cadmium	200.8	200.7	Di(2-Ethylhexyl)Phthalate	525.2
Chromium	200.8	200.7	Diquat	549.2
Copper	200.8	200.7	EDB (Ethylene Dibromide)	504.1
Lead	200.8		Endothall	548.1
Mercury	200.8	245.1	Glyphosate	547
Nickel	200.8	200.7	Oxamyl (Vydate)	531.1
Selenium	200.8		PCB (Polychlorinated Biphenyls)	508.1
Thallium	200.8			
			Disinfection Byproducts	
Inorganics			Bromate	Not Certified
Cyanide	335.4		Chlorite	Not Certified
Fluoride	SM4500 F-C		Haloacetic Acids	552.2
Nitrate	353.2		Total Trihalomethanes	524.2
Nitrite	353.3		Total Organic Carbon	5310C
Volatile Organic Chemicals	524.2		Radionucleides	Not Certified
Pesticides			Dioxin	Not Certified
Alachlor (Lasso)	508.1	525.2		
Atrazine	525.2		Asbestos	Not Certified
Bhc-Gamma (Lindane)	508.1	525.2		
Chlordane	508.1			
Endrin	508.1	525.2		
Heptachlor	508.1	525.2		
Heptachlor Epoxide	508.1	525.2		
Hexachlorobenzene	508.1	525.2		
Hexachlorocyclopentadiene	508.1	525.2		
Methoxychlor	508.1	525.2		
Simazine	525.2			
Toxaphene	508.1			
Herbicides				
2,4,5-TP (Silvex)	515.3			
2,4-D	515.3			
Dalapon	515.3			
Dinoseb	515.3			
Pentachlorophenol	515.3			
Picloram	515.3			